IMPORTANT SAFETY INSTRUCTIONS

- The thermometer is suitable exclusively for measurement of the body temperature!
- The minimum measurement time until the signalling tone (beep) must be maintained without exception!
- Ensure that children do not use this device unsupervised, some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect the thermometer from hitting and striking!
- Avoid ambient temperatures above 50 °C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- The manufacturer recommends verifying the accuracy by an authorised laboratory every 2 years.
- This product is not to be disposed together with other domestic waste after the expiration of service life. To prevent the possible damage to the environment or human health through uncontrolled waste disposal, it is necessary to separate this device from other types of waste and dispose of it in accordance with the recommendations of national standards and rules of disposal, for rational recycling of material resources. This product does not contain harmful substances.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

WARNING: The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only!

TURNING ON THE THERMOMETER

To turn on the thermometer, push the button on the top; a short beep signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, an «-» and a flashing «C» appear at the upper right of the display. The thermometer is now ready for use.

USING THE THERMOMETER

Choose the preferred measuring method. During a measurement, the current temperature is displayed continuously. The «C» symbol flashes during this period of time. If the signalling tone sounds (beep-beep-beep 10 times) and the «C» is no longer flashing, this means that the measured increase in temperature is now less than 0.1 °C in 8 seconds and that the thermometer is ready for reading. To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display. After 3 more seconds, this value will disappear and the thermometer will switch to the normal measurement mode. Then the stored value is erased.

ILLUMINATED DISPLAY / NIGHT GLOW FEATURE

The thermometer will beep and the Night Glow will come ON simultaneously for 5 seconds when you push the on-off button to turn ON the unit, the display will read «188.8» for the first 2 seconds then automatically display the last temperature for the other 3 seconds. When the power button is pressed to Switch-Off the unit, there will be a 4-second additional delay with the Night Glow ON before power off, which allows the user to read the measurement reading easily in the dark. This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

GARANTEE

We grant you a 2 year guarantee after the date of purchase. Any damage caused by improper handling shall not be covered by the guarantee. Battery and packaging are also excluded from the guarantee. All other damage claims excluded. A guarantee claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

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